



**Long Cuff (Elbow Length) Powder Free
Nitrile Examination Gloves (16 Inch & 18 Inch)**

www.kanamlatex.com

Kaltex Nitrile gloves with extra length provides an excellent barrier protection from cross contamination.

SAFETY FEATURES

- ▶ Extra long cuff (400 & 450 mm) covers the gown sleeves, and protects the users forearm against fluid contamination
- ▶ Latex free and Powder free Nitrile examination gloves
- ▶ The micro rough texture of glove surface gives good grip
- ▶ Ambidextrous for quick and easy donning
- ▶ Extra length protects against AIDS/HIV infections
- ▶ Manufactured in a clean environment to comply with EN 556 & ISO 13485 requirements
- ▶ Ideal for pharmaceutical use
- ▶ Beaded cuff for easy donning

STANDARDS

- ▶ ASTM D 6319-2019, EN 455-1:2020 & EN 455-2:2015
- ▶ Certified to meet PPE Standards EN 374, EN 420 - Cat III
- ▶ CE Certified.

QUALITY CONTROL

- ▶ Factory Standard of AQL 1.5



Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppoottil Buildings, K.K. Road
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614

E-mail: customercare@kanamlatex.com
ho@kanamlatex.com

Product Information

TYPE

Powder Free Nitrile Examination Gloves.

PRIMARY MATERIAL

Synthetic Nitrile Rubber Latex

DESIGN

Ambidextrous, Beaded cuff, Micro rough textured.

POWDER FREE RESIDUE

Powder Free: <2mg/glove



PHYSICAL DIMENSIONS

Size	Length (mm) Min	Palm Width (mm) min	Thickness (mm) Min		
			Cuff	Palm	Finger
6 ½	400	80 ± 10	-	0.05	0.05
7 ½		95 ± 10			
8 ½	450	110 ± 10			

PHYSICAL PROPERTIES

Characteristics	Before Ageing	After Ageing 70±2°C for 166±2hrs
Tensile Strength (Mpa) min.	14	14
Ultimate Elongation (%) min.	500	400
Force at break (N) Min.	6.0	6.0

PRE-SHIPMENT QUALITY INSPECTION

Sampling Procedure: ISO 2859 Part I

Sampling Plan: Single sampling, Normal Inspection

Characteristics	Inspection Level	AQL
Freedom from Holes	G-1	1.5
Visual attributes – Major Defects	G-1	1.5
Visual attributes – Minor Defects	G-1	2.5
Physical Dimension	S-2 / N=13	4.0
Physical Properties	S-2	4.0
Force at break	N=13	N/A

PACKAGING - 400 mm

Material	Quantity	Dimension (cm)	
		Sterile	Non-Sterile
Paper Wallet	One Pair	25(L) x 11(W)	N/A
Pouch	One Pair	32(L) x 15(W)	32(L) x 15 (W)
Inner Box	50 Pairs	32(L) x 15(W) x 23(H)	32(L)x15(W)x19.5(H)
Master Carton	250 Pairs	76(L) x 32.5(W) x 24(H)	77(L)x33(W)x21.5(H)

PACKAGING - 450 mm

Material	Quantity	Dimension (cm)	
		Sterile	Non-Sterile
Paper Wallet	One Pair	29(L) x 11(W)	N/A
Pouch	One Pair	36(L) x 15(W)	35(L) x 15 (W)
Inner Box	25/50 Pairs	36(L) x 15(W) x 11.5(H)	35(L)x15(W)x19.5(H)
Master Carton	100/250 Pairs	37(L)x 31(W)x24(H)	77(L)x36(W)x21.5(H)



Long Cuff (Elbow Length) Powder Free
Nitrile Examination Gloves (16 Inch & 18 Inch)