



Long Cuff Powder Free Nitrile Examination Gloves (290mm)

Kaltex Nitrile gloves with extra length provides an excellent barrier protection from cross contamination.

SAFETY FEATURES

- ▶ Long cuff (290 mm) covers the gown sleeves
- ▶ Latex free and Powder free Nitrile examination gloves
- ▶ The micro rough texture of glove surface gives good grip
- ▶ Ambidextrous for quick and easy donning
- ▶ Extra length for extra protection
- ▶ Manufactured in a clean environment to comply with EN 556 & ISO 13485 requirements
- ▶ Ideal for pharmaceutical use
- ▶ Beaded cuff for easy donning

STANDARDS

- ▶ ASTM D 6319-2010, EN 455 -1: 2000 & EN 455-2: 2015
- ▶ Certified to meet PPE Standards EN 374, EN 420 - Cat III
- ▶ CE Certified.

QUALITY CONTROL

- ▶ Factory Standard of AQL 1.5



Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

Oppoottil Buildings, K.K. Road,
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614
E-mail: customercare@kanamlatex.com
kanamlatex@gmail.com

www.kanamlatex.com

POWDERFREE STERILE & NON STERILE NITRILE EXAMINATION GLOVES

Product Information

TYPE

Powder Free Nitrile Examination Gloves.

PRIMARY MATERIAL

Nitrile.

DESIGN

Ambidextrous, Beaded cuff, Micro rough textured.

POWDER FREE RESIDUE

Powder Free: <2mg/glove

PHYSICAL DIMENSIONS

Size	Length (mm) Min	Palm Width (mm) min	Thickness (mm) Min		
			Cuff	Palm	Finger
6 ½ (S)	290	80 ± 10	-	0.05	0.05
7 ½ (M)		95 ± 10			
8 ½ (L)		110 ± 10			

PHYSICAL PROPERTIES

Characteristics	Before Ageing	After Ageing 70±2°C for 166±2hrs
Tensile Strength (Mpa) min.	14	14
Ultimate Elongation (%) min.	500	400
Force at break (N) Min.	6.0 (min.)	6.0 (min.)

PRE-SHIPMENT QUALITY INSPECTION

Sampling Procedure: ISO 2859 Part I

Sampling Plan: Single sampling, Normal Inspection

Characteristics	Inspection Level	AQL
Freedom from Holes	G-1	1.5
Visual attributes – Major Defects	G-1	1.5
Visual attributes – Minor Defects	G-1	2.5
Physical Dimension	S-2 / N=13	4.0
Physical Properties	S-2	4.0
Force at break	N=13	N/A

PACKAGING - 290mm Non Sterile

Material	Quantity	Dimensions (cm)
Inner Box	100 Pcs.	26(L) x 12.5(W) x 7.2(H)
Master Carton	1000 Pcs.	37(L) x 26(W) x 27(H)

PACKAGING - 290mm Sterile

Material	Quantity	Dimensions (cm)
Paper Wallet	1 Pair	19(L) x 10(W)
Pouch	1 Pair	24(L) x 13(W)
Inner Box	50 Pairs	24(L) x 13(W) x 18(H)
Master Carton	300 Pairs	40(L) x 25(W) x 37(H)

