



# PRODUCT SPECIFICATION

BRAND : KALTEX

(400 mm)



## Product : Powder Free Non - Sterile Synthetic Nitrile Examination Gloves

### PHYSICAL DESCRIPTION

Material	Synthetic Nitrile Rubber Latex
Colour	Blue
Shape	Ambidextrous
Cuff	Beaded
External Surface	Chlorinated
Internal Surface	Polymer coated



### PHYSICAL DIMENSION


Size	Minimum Length (mm)	Palm Width (mm)		Weight per glove (g)	Minimum Thickness (mm)		
	KANAM	ASTM D 6319	KANAM		Location	ASTM D 6319	KANAM
Small	400	80 ± 10	84 - 85	12.5 - 13.0	Palm	0.05	0.18
Medium	400	95 ± 10	96 - 97	13.0 - 13.5			
Large	400	110 ± 10	106 - 107	13.5 - 14.0	Finger	0.05	0.19
X-Large	400	120 ± 10	112 - 113	14.0 - 14.5			

### PHYSICAL PROPERTIES

Characteristics	Before Ageing		After Ageing (70 ± 2° C for 166 ± 2 h)	
	ASTM D 6319	KANAM	ASTM D 6319	KANAM
Tensile Strength (Mpa)	14 min	14 min	14 min	14 min
Ultimate Elongation (%)	500 min	500 min	400 min	400 min

<b>Freedom from holes</b>	Inspection Level 1 AQL 1.5	<b>Classification</b>	Class I
<b>Visual-Major Defects</b>	Inspection Level 1 AQL 1.5	<b>Compliances</b>	ASTM D 6319
<b>Visual-Minor Defects</b>	Inspection Level 1 AQL 2.5		ISO 13485, ISO 9001, ISO 14001
<b>Powder Free Residue</b>	<2 mg / glove	<b>Shelf Life Period</b>	3 years
<b>Packaging</b>	10 boxes x 50 pcs / 500 pcs per carton		

**Kanam Latex Industries Private Limited**  
 Site1 : 3/13 F, West Peruvilai Road, Pallavilai - 629 003  
 Site 2 : 12/67C, Ananthanadarkudy - 629 201,  
 Nagercoil, Kanyakumari District,  
 Tamil Nadu, India



Issue # 01, Rev#01  
 Date : 12.01.2021