

Surgicare[®] Aloe Vera

Powderfree Soft Nitrile Examination Gloves

www.kanamlatex.com

Specially designed latex free and powder free non sterile Nitrile Examination gloves, to suit every Healthcare Professional. **The inner coating of Active Aloe Vera moisturizes and cools the user's skin preventing dryness and itching.** It also provides excellent barrier protection for the healthcare worker. Perfect for dentists who have to wear gloves constantly.

SPECIAL FEATURES

- ▶ Inner coating of Active Aloe Vera moisturizes and cools the user's skin preventing dryness and itching
- ▶ Latex free and Powder free
- ▶ Ambidextrous for quick and easy donning
- ▶ Finger textured for additional grip on instruments
- ▶ Weight of 3.5 gms.
- ▶ ASTM : D-6319-2019
- ▶ Device classification - Class 1

Nitrile Exam gloves is the perfect choice for:

- ▶ Dental and ICU Workers
- ▶ Can be used in the food industry and hospitality sector.

QUALITY STANDARDS

- ▶ Manufactured in accordance to ISO 9001:2015 and ISO 13485: 2016 Quality System Standards.

QUALITY CONTROL

- ▶ AQL - 1.5

Moisturizes
& Cools
the Skin



SOFT NITRILE MEDICAL EXAMINATION GLOVES

Kanam
KANAM LATEX
INDUSTRIES PVT. LTD.

The trusted name the world over

Ooppottil Buildings, K.K. Road
Kottayam-686 001, Kerala, India
Tel: 0481-2560108, 2563513
Fax: 0481-2563614

E-mail: customercare@kanamlatex.com
ho@kanamlatex.com

AN ISO 9001:2015, ISO 13485:2016 & ISO 14001:2015 CERTIFIED COMPANY

Product Information

TYPE

Powder Free Nitrile Examination Gloves.

PRIMARY MATERIAL

Nitrile

DESIGN

Ambidextrous, Beaded cuff and textured.

POWDER FREE RESIDUE

<2mg/glove

SHELF LIFE

3 years in original packaging if stored as per storage instructions.

PHYSICAL DIMENSIONS

Size	Length (mm)		Palm Width (mm)		Average wt. per glove (gm)	Minimum Thickness (mm)		
	ASTMD	KANAM	ASTM D	KANAM		Location	ASTM D6319	KANAM
X-Small	220	230	70±10	76-79	3.0-3.2	Palm	0.05	0.07
Small	220	230	80±10	86-89	3.3-3.5			
Medium	230	230	95±10	98-99	3.6-3.8			
Large	230	230	110±10	107-109	3.9-4.1	Finger	0.05	0.10
X-Large	230	235	120±10	114-119	4.2-4.5			

PHYSICAL PROPERTIES

Characteristics	Before Ageing		After Ageing 70±2°C for 166±2hrs	
	ASTM D 6319	KANAM	ASTM D 6319	KANAM
Tensile Strength (Mpa)	14 min	14 min	14 min	14 min
Ultimate Elongation (%)	500 min	500 min	400 min	400 min

PRE-SHIPMENT QUALITY INSPECTION

Sampling Procedure: ISO 2859 Part I

Sampling Plan: Single sampling, Normal Inspection

Characteristics	Inspection Level	Method	AQL
Freedom from Holes	G-1	ASTM D 5151:2011	1.5
Visual attributes – Major Defects	G-1	KL/SP/903	1.5
Visual attributes – Minor Defects	G-1	KL/SP/903	2.5
Physical Dimension	S-2 / N=13	ASTM D 6319:2019	4.0
Physical Properties	S-2	ASTM D 6319:2019	4.0
Extractable Protein Content	N=3	ASTM D 5712:2010	N/A
Powder Amount	N=2	ASTM D 6124:011/ISO 21171:2006	N/A

PACKAGING

Material	Quantity	Dimension (cm)
Box	100 Pieces	23.5(L) x 12.5(W) x 5(H)
Master Carton (20 Inner Boxes)	2000 Pieces	51(L) x 26(W) x 25.5(H)

KLIP/G/B/21/Rev#01



Surgicare[®]

Aloe Vera

Powderfree Soft Nitrile Examination Gloves

